




User Guide

Direct Detect™ Spectrometer Software



Notice

The information in this document is subject to change without notice and should not be construed as a commitment by EMD Millipore Corporation ("Millipore") or an affiliate. Neither EMD Millipore Corporation nor any of its affiliates assumes responsibility for any errors that may appear in this document.

The Direct Detect™ system is intended for research use only and is not for use in diagnostic procedures.

The M logo and Direct Detect are trademarks of Merck KGaA, Darmstadt, Germany.

Millipore is a registered trademark of Merck KGaA.

Microsoft, Excel, and Windows are a registered trademarks of Microsoft Corp. in the United States and/or other countries.

© 2012 EMD Millipore Corporation, Billerica, MA, USA. All rights reserved.

00000844, Rev. A, 02/12

Contents

| | | |
|----------|--------------------------------------|-----------|
| 1 | Software Overview | 1 |
| 1.1 | Program Startup and Login | 1 |
| 1.2 | User Management..... | 2 |
| 1.3 | Main Screen Layout..... | 4 |
| 1.4 | Instrument Operation | 6 |
| 2 | Sample Measurement..... | 8 |
| 2.1 | Reference Measurements | 8 |
| 2.2 | Sample Measurement Setup..... | 9 |
| 2.3 | Sample Measurement | 12 |
| 2.4 | Sample Results..... | 14 |
| 2.5 | Spectrum Plot..... | 15 |
| 2.6 | Curves Plot..... | 18 |
| 3 | Post Analysis..... | 20 |
| 3.1 | Card History..... | 20 |
| 3.2 | Card Statistics | 21 |
| 4 | Administrative Functions..... | 23 |
| 4.1 | Drying Calibration..... | 23 |
| 4.2 | Creating a Calibration Method | 24 |
| 5 | Diagnostics..... | 30 |
| 5.1 | Hardware Diagnostics | 30 |
| 5.2 | Instrument Validation | 33 |
| 6 | File Menu | 35 |
| 6.1 | Print Loaded Record | 35 |
| 6.2 | Logout..... | 36 |
| 6.3 | Exit | 36 |
| 7 | Options Menu | 36 |
| 7.1 | Program Settings..... | 36 |
| 7.2 | Clear Results..... | 39 |
| 7.3 | Select Blank File | 39 |
| 7.4 | Plotting..... | 40 |

Contents, continued

| | | |
|-----------|--------------------------------|-----------|
| 8 | Help Menu | 40 |
| 8.1 | User Guide..... | 40 |
| 8.2 | About..... | 40 |
| 9 | Appendix | 41 |
| 9.1 | Software Updates..... | 41 |
| 9.2 | Folders and Files..... | 41 |
| 10 | Technical Service | 42 |
| 11 | Standard Warranty..... | 42 |

1 Software Overview

Direct Detect™ Spectrometer Software is an application program designed for the Windows® 7 operating system and comes pre-installed on the Netbook computer.

1.1 Program Startup and Login

The Direct Detect™ instrument must be connected to the computer and a power source before starting the Direct Detect™ software (refer section 6.5 in the Direct Detect™ Spectrometer User Guide).

To start the software, click on the Direct Detect™ icon  on the desktop. The software can also be accessed from the "Start" menu and taskbar.

WARNING: Do not install third-party software (screen savers, etc.) as this may impact the performance of the Direct Detect™ system.

During startup the title screen will be displayed. The initialization of the connections to the instrument is configured during the display of the title screen.



Figure 1. Direct Detect™ startup screen

After loading the analysis engine you should hear an audible tone from the instrument to indicate that communication has been established.

If a message appears indicating the software is not able to communicate with the instrument, confirm that the instrument is powered up, that the RJ45 cable is connected between the instrument and PC.

After the software initialization is complete the login dialog appears.




Figure 2. Direct Detect™ login screen

Each user can be assigned a unique name and optional password. The measurements performed by each user are stored and recalled based upon the user name. Certain customization settings are also stored on a per user basis.

1.1 Program Startup and Login, continued

Users can be assigned one of two privilege levels; Operator or Administrator. Administrators have full access to all program features while operators are restricted to sample measurement and analysis only.

By default two users are created. The generic user name "Operator" is a pre-defined operator who has basic sample measurement and analysis capabilities. No password is required for the generic operator.

The "Administrator" user name allows you to login into the **User Management** screen (section 1.2). This screen is used to assign new users and set their privilege level and password.

To log into the program, enter or select a user name from the "User Name" field. Only previously defined users can be entered or selected. The privilege level for the selected user is displayed in braces above the field.

If a password has been assigned to the selected user, then the password field becomes active. You must enter the correct password (not case sensitive).

Press enter or select the "OK" button to login in. After a few moments the main program screen should appear.

Select the "Logout" option from the "File" menu (Ctrl-L) to logout of the main screen and return to the login screen.

1.2 User Management

The **User Management** dialog is used to create, update, or delete users of the Direct Detect™ software. The functions in this section are accessible to users with "Administrator" privilege only.




Figure 3. User management dialog

The **User Management** dialog is accessible from the login screen by selecting the "Administrator" user name. The default password for the administrator is "directdetect".

The **User Management** dialog is also accessible from the main screen by selecting "User Management" from the "Options" menu. This option is only available to users with administrator privileges.

1.2 User Management, continued

The **User Management** dialog has a list of currently defined users on the left side of the dialog, and controls to add, update, or delete users to the right. Left click on a user's name in the list and the current settings for that user appear in the controls on the right.

The "Administrator" password can be modified to prevent access to the **User Management** screen from login. To restore the default password either type it as the password for the administrator, or leave the password field blank and the software will ask you if you want to reset the password to the default.

Adding a new user:

To create a new user, enter a unique name in the "User Name" field. If you want to password protect this user enter a password in the "Password" field, otherwise leave the field blank.

- ☒ Operator Select whether you want this user to have administrative rights which grant access to administrative functions, or to limit the user's access to sample acquisition and analysis.
- ☐ Administrator




Select the "Add" button to add this new user to the list of user names.

Modify a user:



To change a user's password or privilege level, first select their user name from the list, and then make the desired changes. When complete select the "Update" button to save the changes.

Deleting a user:



To delete a user first select the user name from the list and then select the "Delete" button. This will delete the user but not the data files or database created by the user during sample acquisition and analysis.



☐ Show Password Normally the user's password is masked by asterisks. Check the "Show Password" option to display the password.

1.3 Main Screen Layout




Figure 4. Main Direct Detect™ screen

The main Direct Detect™ screen is used to perform all sample measurement analyses. The screen is divided into five regions.




Figure 5. Main screen regions

Sample Information:

This section defines the sampling parameters. Fields and options in this section must be defined before sample measurement can commence.

Sample Measurement:

This section starts the sample measurement and provides feedback to the user during the measurement cycle.

1.3 Main Screen Layout, continued

Measurement Results:

This section contains two tabs. The "Last Sample Results" tab displays the results from the last sample card measured, or the results from the sample card that was displayed using the "Load History" button on the "Card History" tab.

Graphics & Analysis:


This section contains three tabs to display either the spectral plot for a specific region or a calibration curve plot, or a statistics grid of past results.

Diagnostics:

This section contains three buttons. The top button displays the program information. The other two buttons display instrument validation and hardware diagnostics respectively. The icon and tooltip indicate the current state of these diagnostics. Selecting either button opens a dialog with more detailed information.



The top half of the screen can be hidden to allow more screen space for the bottom section. To hide the top half of the screen left-click on the center of the divider. The bottom half of the screen will expand to cover the entire program screen. Click the center divider again to return the screen layout to the default.



The bottom half of the screen also contains a horizontal divider which can be dragged left or right to reveal more or less of each region. Place the mouse pointer over the center of the divider. The pointer will change shape. Left-click and hold the mouse down to drag the divider left or right.

1.4 Instrument Operation

The Direct Detect™ spectrometer measures the near-infrared spectrum of the energy that passes through the protein sample. The chemical composition of the sample is defined by the amount of energy absorbed at specific frequencies, or wavenumbers, of energy.

Samples are introduced using the Assay-free sample card. The card contains four polymer membrane sample positions. An alignment notch in one corner of the card prevents incorrect card insertion into the card holder.



Figure 6. Sample card

By default, the bottom position 1 (closest to the notch) is reserved for the "blank" measurement. The "blank" is the buffer solution without a protein sample. One or more of the remaining positions can be used for samples.

Refer to the section 10.2 in the Direct Detect™ Spectrometer User Guide for detailed instructions on loading and installing the Assay-free sample card. Apply exactly 2 μ L of solution to the blank and sample positions.



Figure 7. Sample card spotting

Insert the card vertically into the slot in the top of the card holder (Figure 8). The card notch should face towards the center of the instrument.

1.4 Instrument Operation, continued



Figure 8. Inserting sample card

When the card is inserted properly the instrument will respond by jogging the card holder up and down briefly. The bottom green LED on the card holder module will light and the system will produce a brief tone. If this does not occur, check the card orientation. The software will not allow sample measurement if the card is not inserted properly.

The top of the card holder module contains two LED lights located in the bottom right corner (as viewed from the top).




Figure 9. Top view card holder module

The bottom LED (green) is lit whenever a sample card is inserted properly.

The top LED (yellow) is lit whenever the drying system heater and fan are active.

Periodically the software will perform a reference measurement which must be performed without a sample card inserted in the card holder. If a card is found in the card holder the software will prompt you to remove it.

Important: A seven minute warm-up cycle is required after initial power-up and software activation. The software will wait for seven minutes before allowing a sample measurement to begin.

2 Sample Measurement

2.1 Reference Measurements

After login to the main screen the software will automatically attempt to perform a reference measurement. This is an internal instrument measurement which subtracts the effects of atmospheric conditions, internal components, etc., from the resultant sample spectrum. A sample card should not be inserted into the card holder while this measurement is in progress.

If there is a sample card in the card holder when an automatic reference measurement begins, a popup dialog will appear requesting that you cancel your measurement so that the reference measurement can begin.




Figure 10. Cancel Card Measurement warning

During the automatic reference measurement a progress dialog will appear. Program access is restricted while the reference measurement is in progress. Do not move, handle, or jostle the instrument while the measurement is in progress.




Figure 11. Reference measurement progress dialog

Once the reference measurement is complete normal program operation can continue. The reference measurement will occur again when the interval defined in the **Program Settings** screen (section 7.1) has expired. The factory-set default setting is 15 minutes.

2.2 Sample Measurement Setup

To perform a sample measurement, prepare the sample card as directed in section 10.2 of the Direct Detect™ Spectrometer User Guide, then complete the fields and options in the sample information region of the main screen.

NOTE: If you choose not to use one of the three preloaded calibration methods, a calibration method must be created prior to reading samples (see Calibration Method, section 4.2).

| Sample Name | Volume μ L |
|-------------------|----------------|
| [Not Used] | 2 |
| Enter Sample Name | 2 |
| sample 1 | 2 |
| blank 1 | 2 |

Figure 12. Sample information region

Each sample card measured requires at least one "Blank" position to be defined. The blank contains the buffer solution used to prepare the samples, but has no protein in it. A previously measured "Blank" can be reused if desired (see option "Use last blank" below). The blank measurement is subtracted from the sample measurement(s) prior to sample analysis.

User Name:


This field is not editable and contains the logged in user name.

Card Name:

Enter a name or identifier for the current sample card.

Method:

This field is a drop-down list of available calibration methods. By default the files are located in "Libraries" (Libraries\Documents\My Documents\Direct Detect\Methods). Three preloaded calibration method files are available (NIST BSA AM1.q3, NIST BSA AM2.q3, and NIST BSA AM3.q3; refer to section 4.2 for details).

If desired use the  button to select a different folder containing Direct Detect™ calibration method files. A folder selection dialog will open to select a new folder.

Below the method field is a description of the method file. Placing the mouse over the description displays the full text of the description in a popup tooltip.

2.2 Sample Measurement Setup, continued

Use last blank:

With this option selected, the last blank collected by the instrument will be used. This allows you to designate all four positions on the card as samples. This option is not available if a previous blank measurement was not performed.

User Name : Ivona

Card Name : Sample1

Method : NIST BSA AM1.q3

NIST BSA standards 0.25mg/ml to 5 mg/ml


☒ Use last blank

☒ Dry sample card

| | Sample Name | Volume µL |
|---|-------------------|-----------|
| 4 | Enter Sample Name | 2 |
| 3 | Enter Sample Name | 2 |
| 2 | Enter Sample Name | 2 |
| 1 | Enter Blank Name | 2 |

Figure 13. Use last blank option

With this option selected an info icon ⓘ will appear. Click on the icon to view details about the blank that will be used.



If you wish to use a previously collected blank other than the last one collected, go to the "Options" menu and choose the "Select Blank File" option (Figure 14). A dialog containing a list of all blank files measured by the current user is presented, with the most recently measured blank at the top of the list.




Figure 14. Select blank dialog

Select an entry from the list and press the **Select** button. The dialog will close and the "Use Last Blank" option will be automatically selected.


2.2 Sample Measurement Setup, continued

Dry sample card:

The option is checked by default and enables the software to perform the card drying routine prior to card measurement and analysis. You can override the "Dry sample card" action by unchecking this option if the sample card has been previously dried, but it is not recommended.

Defining samples:

The right side of the sample information screen contains fields for defining the individual sample positions for the sample card to be measured.





The screenshot shows a table with two columns: 'Sample Name' and 'Volume μL '. There are four rows, each preceded by a colored circle in a square box. The circles are numbered 1 to 4 from bottom to top. Row 1 (blue circle) has 'blank 1' in the name field and '2' in the volume field. Row 2 (green circle) has 'sample 1' in the name field and '2' in the volume field. Row 3 (green circle) has 'Enter Sample Name' in the name field and '2' in the volume field. Row 4 (white circle) has '[Not Used]' in the name field and '2' in the volume field.

| | Sample Name | Volume μL |
|---|-------------------|----------------------|
| 4 | [Not Used] | 2 |
| 3 | Enter Sample Name | 2 |
| 2 | sample 1 | 2 |
| 1 | blank 1 | 2 |

Figure 15. Sample definition fields

Four colored circles represent the four sample card positions. By default the bottom position is designated as a blank position. The color of the circle indicates the use of the sample position. Click on the sample position to change it to the desired designation (blank, sample, or unused).

 A blue circle indicates a blank (buffer) position.

 A green circle indicates a sample position.

 A white circle is an unused position.

Blank and sample positions require a sample name.

If a required field has not been defined, a red circle with an exclamation mark will appear next to the field that is missing when you attempt to start a measurement cycle.

Card positions are measured in numerical order.

Defining a card blank position:

By default the blank position is defined as position #1 on the sample card. The actual position can be defined as any one of the four card positions.

2.2 Sample Measurement Setup, continued

More than one "Blank" position can be defined if desired. Below is an example illustrating the behavior if more than one blank position has been defined.

Example:



The blank position with the lowest position number will be measured first. In this example position 1 is the first defined blank position. If position 3 was the only defined blank then the software would start with position 3 before returning to the lowest to highest numbering scheme.

Position 2 will have the blank from position 1 subtracted from it.

Position 3 will be measured as a blank.

Position 4 will have the blank from position 3 subtracted from it.

2.3 Sample Measurement

After completing the sample measurement setup (section 2.2), and preparing the sample card, you are ready to begin sample measurement.

Insert the sample card into the slot (refer to section 1.4).

The sample measurement region of the screen (Figure 16) contains a colored text status field with a progress bar beneath.






Figure 16. Sample measurement region

Select the  button to start the measurement cycle. After checking for required parameters the measurement cycle will begin.

The system requires a seven minute warm-up period from the time that the instrument is powered before a measurement cycle can be performed. If the warm-up period has not been completed, the status bar will turn yellow  and count down the number of remaining seconds before the measurement cycle can begin.

NOTE: The instrument is powered through the computer and does not have a separate ON/OFF switch. It should be left on and powered up at all times. The seven minute warm-up is required only after the instrument has been powered down.

The Direct Detect™ instrument will position the sample card at the first defined blank position (unless the option to "Use last blank" has been checked) starting at the bottom of the sample card (card position 1).

2.3 Sample Measurement, continued

If the option "Dry sample card" is checked, the instrument will perform the drying routine. The drying fan will start and the yellow LED on the instrument will light.




Figure 17. Drying cycle in progress

The first sample position will typically take the longest period of time to dry, with subsequent positions taking less time.

The status field will indicate which sample position is being dried. If the drying cycle reaches its defined limit for cycles (section 7.1), a warning message will appear.




Figure 18. Drying cycle warning

For samples with a viscosity higher than water (e.g., when glycerol is present), the drying time may exceed the defined time limit. You can continue the drying cycle by selecting the "OK" button, or abort the measurement cycle using the "Cancel" button.

When the drying cycle has finished (unless it was not performed) the sample measurement cycle will commence. The sample data is acquired and analyzed. The progress bar indicates the acquisition progress.





Figure 19. Sample measurement cycle

You can abort the measurement cycle at any time by selecting the  button. The active measurement will be completed before the cycle is aborted and the data will be saved.

The result for each individual sample position is recorded in the "Last Sample Results" tab (section 2.4) of the **Measurement Results** region of the main screen. It can also be viewed and recalled in the "Card History" tab (section 3.1). The spectrum plot and calibration curves are plotted as well (sections 2.5 and 2.6).

2.4 Sample Results

The results from each measured sample position are displayed on the "Last Sample Results" tab. The tab is divided into two sections.

The top section shows the card information, including the name, measurement date and time, and the computed concentration value (mg/mL) for each sample position (Figure 20). The color box to the left of the concentration matches the color of the plotted spectrum and calibration curve point. The sample position name is to the right of the concentration value. For positions that are defined as blank or not used the concentration value is "----".

The bottom section displays the computed average concentration, standard deviation, and % CV (coefficient of variation) for the measured sample positions on the sample card.




Figure 20. Last sample results tab

To remove one or more of the positions from the statistics, hold the control key and left click the sample position circle. The circle color will change to light red (Figure 21). Click again to return it to the measured color and include the result in the statistics.





Figure 21. Last results statistics

When the measurement cycle is complete, select the  button to create a single page report for this sample card. This is the same as selecting "Print Loaded Record" from the "File" menu (section 6.1).

2.5 Spectrum Plot

In the **Graphics & Analysis** region of the screen, the sample card analysis results are displayed as a plot of the spectral data within previously defined regions.




Figure 22. Spectrum plot

Four previously configured spectral regions can be viewed using the tabs along the bottom of the plot. The "Protein" region is selected by default. Select any of the other tabs to view those regions instead.

By default the chart shows the spectral plots for all sample positions measured on the sample card. The blank position(s) can be displayed as well (section 7.1).

The upper-left corner of the plot displays the card name and starting date and time for the card measurement and analysis.

Each sample position is color-coded and matches the color markers displayed on the "Last Sample Results" tab.

Right clicking the mouse on the spectrum plot displays a popup menu of options for the plot.




Figure 23. Spectrum plot menu

Display Legend:



Selecting this option toggles the display of the chart legend. With the legend visible you can hide individual sample plots by unchecking the box to the left of the sample position name.

2.5 Spectrum Plot, continued

Display Toolbar:



Selecting this option toggles display of a toolbar with chart options. This option is available on the "Options" menu as well. Use the hot key combination "Ctrl-T" to toggle this option.



Saves the chart as a bitmap image file.



Copies an image of the chart to the clipboard.



Prints the chart to the default printer.



Displays a printer preview of the chart.

Enable Zooming:

Selecting this option enables chart zooming. Zooming is performed by dragging a box around the area you wish to zoom into, starting from the top-left corner and ending at the bottom-right corner. This option is available on the "Options" menu as well. Use the hot key combination "Ctrl-Z" to toggle this option.




Figure 24. Chart zooming

To reset the zoom press the "Esc" key.

When the chart has been expanded, horizontal and vertical scroll bars appear and allow navigation around the chart.




Figure 25. Chart expanded

Small buttons on the ends of the scroll bars zoom the chart out in discreet increments.

2.5 Spectrum Plot, continued

Display Chart Breaks:

When enabled, the "Display Chart Breaks" option changes the scale of a full spectrum plot by hiding a defined high absorbance region related to the card material (Figure 26.)




Figure 26. Full spectrum without chart breaks (left) and with chart breaks enabled (right)

Display Analysis Region:

Selecting this option highlights the region of the spectrum that is used to compute the concentration values. This is purely a visual aid when examining the spectrum.




Figure 27. Analysis region

2.6 Curves Plot

The "Curves Plot" tab contains a single plotted point for each concentration value computed on the sample card. The color of each point matches the color indicator displayed on the "Last Sample Results" tab (section 2.4).




Figure 28. Curve plot

This plot is based on the calibration method used during sample analysis.

The equation for the line is displayed in the bottom-left corner of the plot and indicates the ideal concentration values per amide signal. The correlation coefficient (R^2) for this method is displayed in the bottom-right corner of the plot.

Zooming for this chart is accomplished by dragging a box around the area you wish to zoom into, starting from the top-left corner and ending at the bottom-right corner.





Figure 29. Curve plot zoomed

To reset the zoom press the "Esc" key.

When the chart has been expanded, horizontal and vertical scroll bars appear and allow navigation around the chart.

Small buttons on the ends of the scroll bars zoom the chart out in discreet increments.



Right click the chart to view a popup menu of options.

2.6 Curves Plot, continued

Display Toolbar:



Selecting this option displays a toolbar with chart options. It is available on the "Options" menu as well, using the hot key combination "Ctrl-T".



Saves the chart as a bitmap image file.



Copies an image of the chart to the clipboard.



Prints the chart to the default printer.



Displays a printer preview of the chart.

Configure Plot Limits:

This option displays a dialog that allows selection of high and low limits for plotting the concentration values.




Figure 30. Plot limits

The default range is 0–10 mg/mL. Enter a new range and select "OK" to apply.

You can reset the range back to the defaults using the "Reset" button.

NOTE: This limit will change automatically if a calibration above 10 mg/mL is performed, or if the protein concentration detected exceeds 10 mg/mL.

3 Post Analysis

3.1 Card History

The "Card History" tab contains two linked data tables which display the raw results from all previous sample measurements for the current user. Each user has a separately maintained database of results.

| Spot # | Sample Name | Concentration | Volume | |
|--------|-------------|---------------|--------|-----------------------------|
| 2 | s1 | 17.277 | 2 | F:\Documents and Settings\F |
| 3 | s2 | 17.282 | 2 | F:\Documents and Settings\F |
| 4 | s3 | 17.283 | 2 | F:\Documents and Settings\F |

| Date | Card Name | Method | Rea |
|----------------------|-----------|----------------------------------|-----|
| 8/3/2011 3:55:46 PM | Sample1 | BSA Standards (Thermo) 071311.q3 | 1 |
| 8/3/2011 3:53:50 PM | Sample1 | BSA Standards (Thermo) 071311.q3 | 1 |
| 8/3/2011 3:35:09 PM | Sample1 | Amide1.q3 | 1 |
| 8/3/2011 3:20:58 PM | Sample1 | Amide1.q3 | 1 |
| 7/26/2011 4:11:20 PM | Sample1 | Amide1.q3 | 1 |
| 7/26/2011 4:03:54 PM | Sample1 | Amide1.q3 | 1 |
| 7/26/2011 4:00:43 PM | Sample1 | Amide1.q3 | 1 |

Figure 31. Card history tab

The bottom table is the list of cards measured with the most recent measurement at the top of the list.


Left click on any record in the bottom table to view the individual result for each sample position on the card in the top table. Only one set of card results can be displayed in the top table.

If a sample spectra file is missing from the hard drive, the file name entry will be highlighted in red.

| Spot # | Sample Name | Concentration | Volume | |
|--------|-------------|---------------|--------|-----------------------------|
| 2 | sample 1 | 1.919 | 2 | F:\Documents and Settings\F |
| 3 | sample2 | 2.041 | 2 | F:\Documents and Settings\F |
| 4 | sample 3 | 1.994 | 2 | F:\Documents and Settings\F |

Right click on either table to view a popup menu of options. For the top table, the options allow you to toggle which columns of data are visible in the table. This selection is user specific and is stored and recalled for the current user.

The bottom table also contains a set of data display filters in the popup menu.



Starting Now:

Clears the list of records from the bottom table and only new records stored after the current date and time will be displayed.


3.1 Card History, continued

Hide Selected:

Hides one or more individual records from view. Hold the control key while left clicking on records to select individual records, or hold the left shift key to select a range of records. Use the hot key combination of "Ctrl-H" to hide the selected records.


Reset Filters:

Resets the actions of the previous two filters.



Create Report

This creates a PDF format report with one or more records depending upon which records were selected.



Load Record

Selecting this button loads the selected record (one only) into the "Last Sample Result" tab, the "Spectrum Plot" tab, and the "Curves Plot" tab.

3.2 Card Statistics

The "Statistics" tab is a table of results from multiple sample cards. Each row in the table is one sample position result.

| | Date | Card Name | Method | Messages |
|---|----------------------|-----------|----------------------------------|----------|
| 1 | 8/4/2011 10:53:44 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |
| 2 | 8/4/2011 10:53:44 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |
| 3 | 8/4/2011 10:53:44 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |
| 4 | 8/4/2011 10:36:43 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |
| 5 | 8/4/2011 10:35:57 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |
| 6 | 8/4/2011 10:28:49 AM | Sample1 | BSA Standards (Thermo) 071311.q3 | |

Average Conc : 1.366 mg/ml

Std Deviation : 0.97

% CV : 70.7 %

Spectrum PlotCurves PlotStatistics

Figure 32. Statistics tab

Adding records to this table generates the same statistics displayed on the "Last Sample Results" tab, but for multiple cards instead of just one card.

The information in the table can be saved to a Microsoft® Excel® format spreadsheet file.

Adding Records to the Table:

Records are added to the table by using the drag-and-drop method of transferring records from either the bottom card history table (to transfer an entire card's worth of data), or by dragging individual sample position results from the top card history table.

3.2 Card Statistics, continued

| | Date | Card Name | |
|---|----------------------|-----------|----------|
| | 8/4/2011 10:36:43 AM | Sample1 | BSA Stan |
| | 8/4/2011 10:35:57 AM | Sample1 | BSA Stan |
| | 8/4/2011 10:30:05 AM | Sample1 | BSA Stan |
| ▶ | 8/4/2011 10:28:49 AM | Sample1 | BSA Stan |
| | 8/4/2011 10:17:30 AM | Sample1 | BSA Stan |
| | 8/3/2011 4:52:48 PM | Sample1 | BSA Stan |

Figure 33. Adding records to the table

Highlight the cards or samples to transfer using the Ctrl and Shift key combinations (Figure 33). Place the mouse pointer at the bottom edge of the last record selected; the mouse pointer will now change shape to include a dotted box below the cursor. Hold the left mouse button down and drag it onto the Statistics table before releasing the mouse button. The record information will be added to the grid.

Editing Records:

Right mouse click the Statistics table to view a popup menu of options.




Figure 34. Statistics popup menu

Displayed Columns:

Toggles which columns of data are visible in the grid.

Export:

Saves the current table as a Microsoft® Excel® format spreadsheet file. A dialog will open to select the file name and location. Only the visible columns are saved to the spreadsheet.

Remove Selected:

To remove rows from the Statistics table, highlight the desired rows using the Ctrl and Shift key combinations, and then select "Remove Selected". The records will be removed from the table (but not from the database).

Remove All:

Select this option to remove all records from the Statistics table.

4 Administrative Functions

All functions in this section are accessible to users with "Administrator" privilege only.

4.1 Drying Calibration

The **Drying Calibration** function is used to calibrate the drying routine used for sample measurement. The default drying routine is adequate for most instruments. However, if this drying routine does not properly dry the sample card, run the **Drying Calibration** function.

From the "Options" menu, choose Drying Calibration (Figure 35).




Figure 35. Drying calibration option

Insert a new (unused), sample card into the Direct Detect™ instrument and select the "Begin Calibration" button. Wait for the drying calibration to complete before removing the sample card.




Figure 36. Drying calibration dialog

When the drying calibration is complete a message will indicate the new drying threshold value that will be used for the instrument. The default value is 0.3.




Figure 37. Drying calibration complete

4.2 Creating a Calibration Method

4.2.1 Calibration Mode

To create a calibration method you will first need to gather sample card spectra from known concentrations of proteins. To gather these files select "Calibration mode (no analysis)" from the "Options" menu (Figure 38). In "Calibration Mode (no analysis)", the main screen is modified to enable sample acquisition with no analysis. These samples are then used to create new calibration method files.




Figure 38. Calibration mode

The calibration methods drop down list is replaced with the "Sample Path" field. This is a folder where all sample files acquired will be stored. Use the button to the right of the field to select the path to store the files. Sample Path :

Because all files measured in this mode will be stored in the same path, it is recommended that the samples names for each sample position be unique.

Prepare and measure sample cards as you would normally. Samples files will be stored in the selected path, and no analyses will be performed on the samples. The results are not displayed and no data is stored in the database for the current user.

To create a robust method you need a number of different samples covering the entire range of concentration values you expect to measure (e.g., 0.2–5 mg/mL). This includes multiple samples that are of the same (or very close) concentration values. A sample set of 20–30 files is generally acceptable for a good calibration method.

4.2.2 Calibration Method

The **Calibration Method** function is used to create or modify calibration method files that will be used to perform the analysis of the sample cards.

Choose "Calibration Method" from the "Options" menu. This opens a wizard that will walk you through the steps of creating or modifying a calibration method file.

NOTE: Calibration method files have the file name extension "q3".

4.2 Creating a Calibration Method, continued

On the first screen of the wizard you decide whether to create a new method or modify an existing method.







Figure 38. Calibration method first screen

The left side of the screen has options for creating or modifying methods. The center of the screen displays a list of method files that exist in the method path. The right side displays information about the selected method file.


The "Path containing method files" field displays the current path to store calibration method files. By default this is located in "Libraries" (Libraries\Documents\My Documents\Direct Detect\Methods).

-  Select the path button to select a different path to store the calibration method files.
-  Use the reset button to return to the default path.


Creating a new calibration method:

To create a new calibration method file, select the  Create New Method option, and then select the "Next >>" button to continue to the next screen of the wizard.

Modifying an existing calibration method:


To modify an existing method file select the file from the "Existing Method Files" list. Then select the option  Modify Existing Method . The spectra files previously used to create this method must still reside in the same path as they existed when the calibration method was first created.

Deleting an existing calibration method:

To delete a method from the list, first select it and then select the  Delete Method button. This will delete only the calibration method file and not the "Calibration Mode" sample files that were used to create the method.

4.2 Creating a Calibration Method *continued*

On the second screen of the **Calibration Method** wizard you select the previously measured sample spectra files and assign the known concentration values for each sample.

Use the  button to open a standard file selection dialog. Select the sample files you want to include in this calibration method. The file names and path are added to the table of files.





Figure 39. Calibration method sample selection

Keep in mind that in order to create a robust method you need a number of different samples (20–30 files) covering the entire range of concentration values you expect to measure.

Select using the Ctrl or Shift keys one or more files to remove from the list and then select the  button.


To remove all the samples from the table, select the  button.

With sample files selected you must enter the known concentration value for the measured sample in the "Conc" column of the table.

When all samples and concentration values have been entered select the "Next >>" button to continue to the next step of the calibration method wizard.

4.2 Creating a Calibration Method, continued

The "Calibrate Samples" screen of the Calibration Method wizard calculates the ideal fit for the samples you selected.

Select the  button to begin the line fitting routine.

When the calibration is complete the results are listed in the table, as well as being plotted against the ideal line.




Figure 40. Calibration method calibrate screen

In a perfect calibration every point would intersect the line exactly. The correlation coefficient listed about the plot (R^2) would be "1.000".

In a good calibration each sample point would be on or near the ideal line. The plot in Figure 40 has a very good correlation coefficient (0.9913), while the plot in figure 41 has a poor correlation coefficient (0.4404).




Figure 41. Bad correlation plot

Typically there will be one or two outliers or samples that are far from the ideal line. This can be caused by a bad sample measurement or an incorrect "Conc" value entered in the previous wizard screen.

Use the "<< Back" button to return to the previous wizard screen. Remove the bad samples and try to calibrate again.

4.2 Creating a Calibration Method, continued

Above the table of concentration/signal values, you can select one of three possible analysis types.






Figure 42. Analysis type

Each analysis type corresponds to a calibration curve (NIST BSA AM1.q3, NIST BSA AM2.q3, or NIST BSA AM3.q3) appropriate for the type of sample analyzed. The three curves compensate for spectral differences and buffer interference by applying different modes of integration. The graphic shows how the peak region for each mode is integrated.




NIST BSA AM1.q3

Method 1 relies on amide peak intensity relative to the local baseline (between predetermined wavenumbers). This method should be used for the analysis of samples where there is no buffer interference and the amide peak is well defined.



NIST BSA AM2.q3

Method 2 allows concentration determination in the presence of interfering buffers (e.g., Tris). The method uses amide signal intensity at a predetermined wavenumber (local baseline established as in Method 1). This method should be used for analysis of samples in buffer that contains interfering components.



NIST BSA AM3.q3

Method 3 compensates for spectral inconsistency to the right of amide signal. Signal intensity is measured at a predetermined wavenumber and the baseline is 'anchored' at a higher (to the left) wavenumber. This method is recommended for quantification of peptides.

Select the analysis type that best corresponds to the analyzed samples.

When you are satisfied with the results of the calibration select the "NEXT >>" button to continue to the next screen of the calibration method wizard.

4.2 Creating a Calibration Method, continued

The final screen of the **Calibration Method** wizard stores the method file as a file.





Figure 43. Calibration method save screen

A list of existing methods is displayed for reference.


If you are creating a new method you must enter a unique name for the method file (file naming restrictions apply).

If you are modifying an existing method then you cannot modify the method name.

In either case you can add or modify a description for this method. Use multiple lines for your description.

Use the  button to save or overwrite the calibration method file.

Select the  button to exit the **Calibration Method** wizard.

Select the  button to return to the start of the wizard to begin another calibration.

5 Diagnostics

For additional information on diagnosing instrument issues, refer to section 11, Troubleshooting, in the Direct Detect™ Spectrometer User Guide.

5.1 Hardware Diagnostics



The hardware diagnostics icon provides a quick visual indication of the state of the instrument's hardware diagnostics. When no warnings or errors are present the icon is green with a green checkmark.



A yellow hardware diagnostics icon indicates some warning condition is present. Under most conditions this does not prevent the user from performing sample acquisition and analysis. If the warning indicates that the instrument warm-up period is active then measurement will not be allowed until the warm-up is complete.



A red hardware diagnostics icon indicates a critical failure in a hardware component. A source error can usually be corrected by replacing the NIR source element. Other types of hardware errors may require instrument repair.



Hardware Diagnostics
Preventative Maintenance Requested
Click for more information

Placing the mouse cursor over the icon displays a tooltip with simple information about the hardware diagnostics state.

Click on the icon to display the hardware diagnostics dialog which contains more detailed diagnostics information and allows you to create a hardware diagnostics report.

The hardware diagnostics dialog displays a row of icons indicating the state of the various hardware components of the instrument. Each icon indicates whether the component is good, reporting a warning, or reporting an error.

5.1 Hardware Diagnostics, continued




Figure 44. Hardware diagnostics source warning

Figure 44 shows the hardware diagnostics screen when the main icon is displaying a warning. In this example the instrument's NIR source is indicating that it is nearing the end of its operational lifetime. Sample measurements can still be made and acquired however performance may suffer and the values for concentration may become suspect. If you receive this warning the source element should be replaced as soon as possible. Refer to section 12.1 of the Direct Detect™ Spectrometer User Guide for instructions on replacing the NIR source.

In Figure 45 the hardware diagnostics is reporting an error with the source and the status message indicates that the source is not connected. This generally means that source element is broken and needs to be replaced.




Figure 45. Hardware diagnostics source error

Figure 46 indicates a warning involving the interferometer components. In this example the status message is indicating that the internal humidity of the instrument is too high. This typically means that the internal desiccant pack needs to be replaced.

5.1 Hardware Diagnostics, continued






Figure 46. Hardware diagnostics interferometer warning

Select the save button  to save a hardware diagnostics report file to the hard drive.

Select the print button  to print the report (Figure 47) to an installed printer.




Figure 47. Hardware diagnostics report

5.2 Instrument Validation



The instrument validation icon provides a quick visual indication of the state of the instrument's operational validation status. When no warnings or errors are present the icon is green with a green checkmark.



A yellow validation icon indicates one of the validation tests has expired. Sample measurement and analysis is still possible when a validation warning is present, however your laboratory practices may require that a new test be performed before continuing to measure samples.



A red validation icon indicates one of the validation tests failed its last measurement. Sample acquisition and analysis is not allowed if a failure has been reported. A failure can indicate the need for instrument service or repair.



Placing the mouse cursor over the icon displays a tooltip with simple information about the validation state.

Click on the icon to display the instrument validation dialog which contains more detailed information.

Figure 48 shows the **Instrument Validation** dialog. The dialog displays identical controls for two different types of tests. Both tests display the status of the tests. Typically this is the amount of time until the test expires, or the amount of time since the test expired.





Figure 48. Instrument validation dialog

5.2 Instrument Validation Diagnostics, continued

The **Performance Qualification (PQ)** test is a test that is typically performed daily. This test compares the performance of the instrument from day to day. A failure of this test may indicate a failure or performance drop in one of the system components (e.g., the source may be degrading).

The **Operational Performance (OQ)** tests the basic performance characteristics of the instrument. This test is typically performed annually. A failure of this test can indicate a hardware misalignment or failure.

A user with administrative privileges can perform these tests using the  buttons. Select either or both of the tests to perform. If you select the other test while one test is active, the software will start the second test as soon as the first test is complete. You can choose to abort the test before it begins.

Important: Do not insert a sample card while either of the validation tests is running.

Figure 49 shows the PQ test running. The status field turns yellow and displays which portion of the test is currently in progress. A progress bar indicates the percentage to completion for the entire test.






Figure 49. PQ test active

To abort the test in progress, select the  button. You will have to wait until the current segment of the test has completed. Generally it is not advisable to abort a test in progress since this will result in a test failure condition which will prevent sample measurement and acquisition.

If you selected a second test to start while the first test is active, and then abort the second test before it has begun, then the test will not be run and the status will remain the same.

5.2 Instrument Validation, continued

You can select the  button to display the last PQ or OQ test report. Reports are stored in the program installation directory under the path "Validation\Reports". This button will display only the most recent report (Figure 50).

| OVP - PQ Test Protocol | | | |
|--|--|---|---|
| Company: | | MILLIPORE | |
| Operator: | | Default | |
| Instrument Type: | | Alpha Sample Compartment RT-DLaTGS ZnSe | |
| Optics Configuration: | | Fibre 1 with: MIR, KBr, RT-DLaTGS (Internal Pos.1) | |
| Accessory: | | Millipore Direct-Detect | |
| Instrument Serial Number: | | M3 | |
| Instrument Firmware Version: | | 1.320_MDD_01 Oct 27 2011 | |
| OPUSDB Version: | | OPUS 7.0 Build: 7, 0, 122 / DB: 7.0.105.1106 for OPUS 7.0 | |
| Overall Test Result | | PASSED | |
| Test expires: | | 12/25/2011, 3:12:41 PM (GMT-5) | |
| Test Date/Time: | | 12/15/2011, 3:12:41 PM (GMT-5) | |
| Test Spectra Path: | | C:\Direct Detect\Validation\Data\20111215\151241 | |
| Date of last PQ Reference Measurement | | | |
| Comment: | | | |
| Signal to Noise Test | | | |
| Minimum Signal to Noise Limit(area 1): | | 1000 | Signal to Noise measured: 3458 |
| 100% Line Test | | | |
| Maximum 100% Line Deviation: | | 2.0 | Measured 100% Line Deviation: 0.29 |
| Wavenumber Accuracy Test - Polystyrene (35.5 Deg. C. Peak is T. corr.) | | | |
| Sample Material: | | Polystyrene | |
| Specified Peak: | | 1601.45 cm ⁻¹ | Maximum Deviation: 0.50 cm ⁻¹ |
| Measured Peak: | | 1601.15 cm ⁻¹ | |
| Corrected Peak: | | 1601.25 cm ⁻¹ | Measured Deviation: 0.20 cm ⁻¹ |
| Overall Test Result = PASSED | | | |
| Date and Signature | | Date and Signature | |

Figure 50. Validation PQ test report

6 File Menu

This section explains the options on the "File" menu of the main program screen.

6.1 Print Loaded Record



Figure 51. File menu

This option prints the single sample card measurement report for the currently loaded sample results. This is either the last measurement and analysis performed, or the result loaded from the database using the "Load Record" button (section 3.1).

Use the hot key combination Ctrl-P to select this option, or use the "Create Report" button on the "Last Sample Results" tab (section 2.4).

6.1 Print Loaded Record, continued

Selecting the "Print Loaded Record" option displays a report options dialog.




Figure 52. Report options dialog

You can print the report to a printer or file (as a PDF report), or "print preview" the report. You can also choose which plot to include with the report; either the spectrum plot or the calibration curve plot.

Select the  button to print to the selected location.

6.2 Logout

Select this option (Ctrl-L) to exit the main program screen and return to the login screen. This is recommended whenever the user plans to leave the instrument unattended for long periods, or to switch users.

6.3 Exit

Select the exit option to close down the program without returning to the login screen.

7 Options Menu

7.1 Program Settings

The options menu is divided into two sections. The bottom section is available only to users with administrator privileges. Refer to section 1.2 for information about administrative functions.




Figure 53. Options menu

7.1 Program Settings, continued

The **Program Settings** option displays a dialog with two tabs; one for user specific settings, and one for global (all users) settings (Figure 54). Each user has settings that can be customized for them.




Figure 54. Program settings – User

Number of measurements to average:

This is the number of consecutive measurements that will be performed for each sample position measured. By default only one sample acquisition is performed. Increasing this number will result in multiple acquisitions which will be averaged before the analysis is performed.

Display blank measurements in history:

Normally the blank positions are not displayed in the top section of the "Card History" tab. Select this option to display the blank positions as well as the sample positions. This will also result in the blank measurement being plotted to the "Spectrum Plot" tab.

Reset "Do not show ..." messages.

Some message boxes will have an option to not show the message again. Select this button to reset these messages so that they will be shown.

Only users with administrative privileges can modify the settings in the **Global Settings** tab of the **Program Options** dialog (Figure 55).




Figure 55. Global settings Options 1

7.1 Program Settings, continued


The global settings are divided into two tabs labeled "Options 1" and "Options 2".


Root Storage Path for Files:

This is the root path for storing all sample files measured during normal acquisition and analysis.

The default path is located in the "Libraries" folder.

"Libraries\Documents\My Documents\Direct Detect\"

Use the selection button  to select a different root path if desired. The actual path for stored sample files depends upon the user and date and time of the measurement (section 9.2).

Use the reset button  if you made changes to the path and want to return to the default path.

Drying Cycle:

The drying cycle section contains the parameters for performing the drying cycle part of the measurement.

Threshold: This is the threshold for the drying routine. Normally this is set by using the "Drying Calibration" function of the "Options" menu.

Cycles: The number of consecutive cycles that must be dry before the routine succeeds.

Attempts: The number of cycles that will be performed before the user is warned that the drying routine is not complete. You can increase this value if you routinely require more than one cycle to complete the drying.

Reference Measurement:

The system performs a reference measurement when the main screen is first displayed, and then again at periodic intervals.

Interval between reference measurements: This is the interval (in minutes) between each successive reference measurement.

Measure: This button allows you to perform a reference measurement immediately. This is useful if you change the interval and want to start the measurements immediately.

Measurement:

Number of Scans: Increase the number of scans to increase the signal-to-noise ratio of samples. This will result in longer acquisition times. It is not recommended that you reduce the number of scans below the default value of 6.

7.1 Program Settings, continued

The "Options 2" tab (Figure 56) contains a table which allows you to customize the tabs displayed on the "Spectrum Plot" tab of the main screen.




Figure 56. Global settings Options 2

A default set of tabs and ranges are provided. Use the "Reset to defaults" button to reset this table to the defaults at any time.

To disable any of the defined tabs uncheck the box in the "Enabled" column.

To delete a row permanently, select the desired row(s) and press the keyboard "Delete" key.

Enabled: Check this column to display this tab.

Tab Name: Name that appears on the tab.

Title: Title for the plot.

X-Min: Minimum x-axis range of plot.

X-Max: Maximum x-axis range of plot.

Interval: The increments between intervals along the x-axis.

A maximum of ten tabs can be defined.

7.2 Clear Results

Select this menu option to clear the "Last Sample Results" and the spectrum and curves plots.

7.3 Select Blank File

When selecting the "Use Last Blank" option on the main screen (section 2.2), the last measured blank is used as the blank measurement for all future measurements. Use this option to select a previously measured blank file from a list.

7.3 Select Blank File, continued

A dialog containing a list of all blank files measured by the current user is presented, with the most recently measured blank at the top of the list.





Figure 57. Select blank dialog

Select an entry from the list and press the  button. The dialog will close and the "Use Last Blank" option will be automatically selected.

7.4 Plotting

The options on this menu mimic the options presented when right-clicking on the "Spectrum Plot" tab. Refer to section 2.5 for more information.

8 Help Menu

8.1 User Guides

Select "Software User Guide" from the "Help" menu or use the hot key combination Ctrl+F1 to display the software user guide, or select "Direct Detect User Guide" (hot key combination Ctrl+F2) to display the spectrometer user guide.





Figure 58. Help menu

8.2 About

Select "About" from the "Help" menu to display a dialog containing information about the software and hardware.

This dialog can also be accessed using the logo button  in the upper right corner of the main screen.

9 Appendix

9.1 Software Updates

Direct Detect™ software may be updated periodically. For information on the most up-to-date software, go to www.millipore.com/directdetect. Register your instrument in order to receive notification about relevant software upgrades.

9.2 Folders and Files

Assorted files and folders are created by the software during sample acquisition and analysis. A defined structure must be maintained for the software to be able to function properly. The sample files can be backed up and restored as needed as long as the entire folder/file structure is maintained.

By default all files and folders are stored in the root file storage path defined in the global settings (section 7.1).

Folders:

By default the root file storage path is in the "My Documents" folder in "Libraries" (Libraries/Documents/My Documents/Direct Detect).




Figure 59. File and folder structure

Folders are defined follows:

Direct Detect: Under the root folder "My Documents" is the folder called "Direct Detect". All files and folders created by the software will be created in this folder.

Note: If an error log is generated it will be stored in this folder.

Methods: The methods folder contains the calibration method files used for sample analysis. These files have a "q3" file extension.

Spectra: All sample files are stored in the "Spectra" folder.

Operator: A sub folder is created for each user as they log into the program. The folder name is the same as the user name.

9.2 Folders and Files, continued

20120208: Each sample set is stored under a subfolder which is based upon the date of the sample acquisition (YYYYMMDD). In this example the date of acquisition was February 8, 2012.

132818: The actual sample files for each sample card measured are stored in a time-based folder. The time is the start of the sample measurement cycle. In this example the measurement cycle started at 13:28:18 (HH:MM:SS). The time value is 24 hour format, so 1:28 PM is 13:28.

Files:

DDUserData.mdb: This is the user-specific database file which contains all the records displayed on the "Card History" tab on the main screen.

LastBuffer.0: This is a copy of the last blank measured. This is what is used when the option to "Use Last Blank" is selected.

Sample 1.2.0: The sample and blank file names are divided into three sections separated by decimal points:

{Sample Name}.{Card Position}.{Type}

Sample Name = The sample name defined by the user for this sample position.

Card Position = The sample card position for this sample (1-4).

Type = '0' is the original file and '1' is the file modified during analysis.

For blank measurements the 'Type' is always '1'.

10 Technical Service

For more information, contact the office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

11 Standard Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms (within the "Terms and Conditions of Sale" applicable to your purchase transaction).